

SOS-Extraction set



SOS set

Extraction damaged bone screws

Compact. Universal. Essential.

endocon[©]

SOS set

Universal op removal tool for damaged and broken bone screws.

SOS set complete consisting of:


All tools available for screw size:

1,5/2,0/2,7/3,5-4,0/4,5/5,0-6,5-7,0 mm

- 1 x Concial extraction screw 1,5/2,0 mm diam.
- 1 x Concial extraction screw 2,7/3,5/4,0 mm diam.
- 1 x Concial extraction screw 4,5/6,5/7,0 mm diam.
- 6 x Extraction bolt
- 6 x Hollow reamer
- 6 x Central bolt for hollow reamer
- 6 x Tube for hollow reamer
- 1 x Hook sharp, lenght 155 mm
- 1 x Screw holding forceps, 20,5 cm / 8"
- 1 x Mini-lexer osteotome, 4 mm
- 1 x Mini-lexer osteotome, 6 mm
- 1 x Mini-lexer osteotome, 10 mm
- 1 x T-handle canullated for AO-shaft 160 mm
- 1 x Aluminium tray for SOS set

SOS-Set - For a smooth operation flow

**Compact.
Universal.
Essential.**



Simplify Mobility

Equipment and accessories

Article no.	Description
54-09-25-0000	SOS set incl. tray



Take advantage of our 24h instrument hire service! We provide all equipment for a fee.

For detailed information on all the endocon products, please refer to the instructions for use or ask our medical product advisors.

This document does not replace an instruction in the intended use of the product.

Distribution:

endocon GmbH Germany
Carl-Bosch-Straße 4 • D-69115 Heidelberg
T +49 6221 7575 55 • F +49 6221 7575 45
E info@endocon.de • I www.endocon.de



endocon